

ADMECO-Lux LED 2KFF

- 60'000 lux at a working distance of 1 m
- LED technology guaranteeing high-energy efficiency and light yield
- Excellent depth illumination
- Dimmer with continuous light regulation
- Color temperature at choice in 3 models (4500K, 3900K oder 3400K)
- Color rendering index > 95
- Light completely free of infrared radiation (no heat radiation)
- Lined casing for optimized projector cooling (hand warm)
- Simple control unit on projector

ADMECO-Lux LED 2KFF

- 60'000 Lux in einer Arbeitsdistanz von 1 m
- Hohe Energieeffizienz und Lichtausbeute dank LED-Technologie
- Hervorragende Tiefenausleuchtung
- Stufenlose Dimmbarkeit der Lichtintensität
- Farbtemperatur wählbar in drei Ausführungen (4500K, 3900K oder 3400K)
- Farbwiedergabe-Index von > 95
- Keine Infrarot-Strahlung im Licht (keine Wärmestrahlung)
- Doppelwandiges Gehäuse zur optimierten Kühlung des Leuchtenkopfes (handwarm)
- Einfache Bedieneinheit am Leuchtenkopf



ADMECO-Lux LED 2KFF-D (also available as wall and mobile version)

ADMECO-Lux LED 2KFF-D (auch als Wand- und Mobilversion erhältlich)

Technical Data	
Illumination intensity at 1m in lux / Beleuchtungsstärke in Lux, 1 m Abstand	60'000
Number of light heads for each model / Anzahl Strahler pro Leuchte	2
Electric power in watts / Elektrische Leistung in Watt	48
Weight in kg / Gewicht pro Leuchte in kg	3.2
Dimmer with continuous light regulation / Stufenlose Dimmbarkeit	yes

Right to change specifications reserved